

Supplementary table 1. Primers used to perform PCR, Sanger sequencing and High Resolution Melting.

SNP	PRIMER SEQUENCES	TECHNIQUE
rs2910164 (<i>pre-miR-146a</i>)	Forward: 5'-ATAACTCATGAGTGCCAGGAC-3'	PCR and sequencing
	Reverse: 5'-TCACAGGAACTCACACTCCTT-3'	PCR and sequencing
rs1862391 (<i>PTTG1</i>)	Forward: 5'-TAGAAACAGTGTCTTAGCCTGA-3'	PCR and sequencing
	Reverse: 5'-AAGCTGCATGGACTGAGGT-3'	PCR and sequencing
	Forward: 5'-ACAACCAATGTTGGGTCTAAA-3'	HRM
	Reverse: 5'-AAACGCTGCTGTTTGATCTA-3'	HRM
rs2910202 (<i>PTTG1</i>)	Forward: 5'-CCCAGATGAGAGGACTTTGCTGT-3'	PCR and sequencing
	Reverse: 5'-GGGCATGAATCACCCAGGA-3'	PCR and sequencing
	Forward: 5'-GGTGAGATAAAGCAGGAA-3'	HRM
	Reverse: 5'-AGGTCTACTACAGTCTACTAT-3'	HRM